ORTHO
●●ONE

REF 161 400 | REF 161 400H

KNIE BANDAGE

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QUALITY. MADE FOR YOU.

ILLUSTRATIONEN
FIGURES

REF D161 400


1


2

3

4

REF D161 400H


1

2

3

4

SPRACHEN LANGUAGES

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
Many thanks for placing your trust in this medical device. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

INTENDED PURPOSE

This medical device is a knee support for compression of soft tissues with a local pressure pad around the kneecap consisting of flexible material.

INDICATIONS

Chronic, post-traumatic or post-operative irritations of knee joint soft tissues, recurrent articular effusion, e.g. in the case of gonarthrosis, arthritis, ligament instability, meniscopathy.

CONTRAINDICATIONS

Thrombosis risk, extreme varicosis, peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders of

the region of the body treated, skin diseases in the part of the body treated, latex allergy.

APPLICATION RISKS / IMPORTANT NOTES

This medical device is a prescribable product. Discuss the use and duration with your treating physician.

The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- remove it during longer periods of rest, e.g. when sleeping
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs

- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use. This medical aid is intended for treating one patient only
- whilst wearing the support, please neither use any creams nor ointments on or around it as they can destroy the material

FITTING INSTRUCTIONS

161 400 Knee Support without a Grip Top

- 1** Grasp the upper edge of the support.
- 2** Pull it over your foot up to the knee with the pad facing forwards.
- 3** The pad should be placed comfortably over the kneecap.

REMOVAL

- 4** To remove the support, please grasp the lower edge at the side on the splints and pull it downwards over the foot.

161 400H – Knee Support with a Grip Top

- 1** Fold the grip top outwards, grasp the upper edge of the support.
- 2** Then pull it over your foot up to the knee with the pad facing forwards.
- 3** The pad should be placed comfortably over the kneecap. Then, flap the grip top back towards the leg so that the silicone nodules are directly on the skin.

REMOVAL

- 4** To remove the support, please fold the grip top outwards, grasp the lower edge of the support at the side on the splints and pull it downwards over the foot.

EN

MATERIAL COMPOSITION

Elastodiene/latex (LA), lyocell (CLY), polyamide (PA), elastane (EL), silicone (SI)

The sewn in textile label provides the precise material composition.

The product contains latex and can trigger allergic reactions.

CLEANING INFORMATION

 Delicate wash 30 °C  Do not bleach

 Do not dry in a tumble dryer  Do not iron

 Do not clean chemically

Do not use fabric conditioner. Stretch back into shape and dry in the fresh air.

GUARANTEE

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions

for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

USEFUL LIFE/LIFETIME OF THE PRODUCT

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

DUTY OF NOTIFICATION

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find the manufacturer's contact data in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link:
www.bort.com/md-eu-contact.

DISPOSAL


Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

DECLARATION OF CONFORMITY

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: www.ortho-one.net/konformitaet

Status: 10.2023

MD Medical device |  Single patient – multiple use

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BORT GmbH