ORTHO  
●●ONE

REF 161 720


# ELLENBOGEN BANDAGE PLUS

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
QUALITY. MADE FOR YOU.

ILLUSTRATIONEN  
FIGURES

1


2

## SPRACHEN LANGUAGES

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**EN** english    Instructions for use ..... 08

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
Many thanks for placing your trust in this medical device. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

### INTENDED PURPOSE

This medical device is a circumferential compression elbow compression support with elastic material fabric pads and elastic profile pads.

### INDICATIONS

Chronic, post-traumatic or post-operative soft tissue irritation in the elbow joint area, epicondylitis, tendomyopathy.

### CONTRAINDICATIONS

Peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders with the region of the body treated, skin diseases in the part of the body treated, latex allergy.

### APPLICATION RISKS/IMPORTANT NOTES

This medical device is a prescribable product. Discuss the use and duration with your treating physician.

The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- remove it during longer periods of rest, e.g. when sleeping
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds

- do not use in case of intolerance of one of the materials used
- no re-use – this medical aid is intended for treating one patient
- whilst wearing the support, please neither use any creams nor ointments on or around it as they can destroy the material

## FITTING INSTRUCTIONS

Place the support on the elbow affected and position the knitted comfort zone in the antecubital fossa at the same time. The support is correctly positioned if you can feel the outer and inner condyle through the hollows in the silicone insert.

**1** Stabilo strap with the easyClick system for individual dosage of the pad pressure: **1** and **3** attachment module, **2** supporting module, **4** adjustment module.

Adapt the pad pressure dosage individually by moving the adjustment module **4**.

**2** When you receive your support, the attachment Stabilo strap module **1** in the supporting module **2** is attached on the support.

Guide the Stabilo strap circumferentially around the lower arm under slight tension.

Close the strap by attaching the free attachment module **3** to the adjustment module **4**. Pay attention to comfortable fitting of the Stabilo strap, neither too loose nor too tight.

If the Stabilo strap with the easyClick system is not needed, the supporting module **2** can be taken from the velour bag.

## REMOVAL

After opening the Stabilo strap, take hold of the lower edge of the support and slip the support over the lower arm in the direction of the hand.

## MATERIAL COMPOSITION

Polyamide (PA), elastodiene/latex (LA), viscose (CV), elastane (EL)

The sewn in textile label provides the precise material composition.

The product contains latex and can trigger allergic reactions.

## CLEANING INFORMATION

 Delicate wash 30 °C  Do not bleach

 Do not dry in a tumble dryer  Do not iron

 Do not clean chemically

Do not use fabric conditioner. Stretch back into shape and dry in the fresh air.

## GUARANTEE

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions

for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

## USEFUL LIFE/LIFETIME OF THE PRODUCT

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

## DUTY OF NOTIFICATION

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find the manufacturer's contact data in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: [www.bort.com/md-eu-contact](http://www.bort.com/md-eu-contact).

## DISPOSAL


Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

## DECLARATION OF CONFORMITY

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: [www.ortho-one.net/konformitaet](http://www.ortho-one.net/konformitaet)

Status: 02.2022

Medical device |  Single patient – multiple use

ORTHO  
ONE

PDF: [ga.ortho-one.net](http://ga.ortho-one.net)

**BORT GmbH**  
Am Schweizerbach 1  
D-71384 Weinstadt  
[www.bort.com](http://www.bort.com)

QUALITY. MADE FOR YOU.

# ELLENBOGEN BANDAGE PLUS

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